

FILAMATIC NEWS

NATIONAL INSTRUMENT COMPANY, INC. • 4119 Fordleigh Road • Baltimore, MD 21215 • Tel. (410) 764-0900 • Fax (410) 764-7719

Internet: <http://www.filamatic.com> • Email: info@filamatic.com

*Manufacturers of Liquid
Filling, Capping and
Packaging Systems*



Contact: Sales Department,
for more information

MULTI-FUNCTION ASEPTIC LIQUID FILLING MONOBLOC SYSTEM **Fill, Purge & Stopper 2mL to 500mL Vials, Bottles** **and Syringe Cartridges on one machine**

The **FILAMATIC**[®] RFS-800 is a multi function, continuous motion, servo controlled filling and stoppering monobloc designed for aseptic liquid filling. Numerous filling methods are available, including frictionless rolling diaphragm pumps for zero-particulate operation. With fill accuracies in the range of $\pm 0.5\%$ or better, it is ideal for applications where both accuracy and sterility are primary concerns.

Incorporation of enhanced self-diagnostics minimizes downtime. Containers are continuously tracked through all machine functions. Fallen bottle reject and sample collection trays are incorporated into the design of this machine. Only precisely filled and stoppered containers are released to the discharge conveyor.

The RFS-800 features a state-of-the-art proprietary Quality Assurance System for Statistical Process Control. A smart system capable of self-correcting fill weights minimizes operator involvement, substantially reducing risk of contamination. Net weight protocols feature automatic feedback that enable individual servo controlled tilling drives to self adjust to pre-established upper and lower control limits, thereby ensuring automatic fill weight compliance.