

FILAMATIC NEWS

NATIONAL INSTRUMENT COMPANY, INC. • 4119 Fordleigh Road • Baltimore, MD 21215 • Tel. (410) 764-0900 • Fax (410) 764-7719

Internet: <http://www.filamatic.com> • Email: info@filamatic.com

Contact: Sales Department,
for more information

*Manufacturers of Liquid
Filling, Capping and
Packaging Systems*



FILAMATIC® Model CE-AB-5 Liquid Filling Machine Now Available With CE Mark

National Instrument Company, Inc., a leading manufacturer of liquid filling, stoppering and capping equipment for the pharmaceutical, cosmetic and chemical specialty industries, has redesigned its basic benchtop Model AB-5 Filling Machine to meet CE standards for the European market.

By adhering to strict European guidelines, National Instrument Company modified **FILAMATIC®** Model AB-5 Filling Machine to Model CE-AB-5 Filling Machine. The IEC (International Electro-Technical Committee) forms the international standards in Europe and provides certain laws for manufacturers to fulfill in order to market their product in Europe. In order for a machine to be endorsed with the CE mark, it has to be designed to meet these specific IEC requirements, including certain safety and electrical compatibility rules. In addition, the machine has to be tested to make sure it passes inspection. To comply with IEC standards, the Model AB machine has undergone some necessary modifications to qualify for the CE mark.

An independent agency was hired to test and evaluate the new CE-AB's

performance. This agency reported, "The National Instrument Model AB-5 Liquid Filling Machine was determined to be compliant with the IEC directives," thus confirming it was ready to sell in the European market.

The new **FILAMATIC®** Model CE-AB-5, a light-duty, bench-top liquid filling machine used to dispense pharmaceutical solutions, cosmetic products and food additives. Depending on which Filling Unit is used, the machine will dispense volumes of free flowing water-thin liquids from .01cc to 130cc per cycle, up to 35 times a minute. The amount of product dispensed from each Filling Unit is consistent and exact, filling volumes accurately to a hundredth of a gram. Operation is simple. The operator can hold a container under the nozzle and press a foot switch to actuate the filler or depress a button on the nozzle bracket to cycle the machine.

National Instrument Company guarantees customer satisfaction and is proud to bear its first CE mark by introducing the new **FILAMATIC®** Model CE-AB-5 Filling Machine to users around the world.