



NATIONAL INSTRUMENT, LLC
Manufacturers of Liquid Filling, Capping and Packaging Systems

FILAMATIC
SPECIALIST IN LIQUID FILLING SYSTEMS

RELEASE AT WILL
Contact: Deven Wlotko
1.410.951.2067 or 1.866.258.1914
Email: info@filamatic.com

February 2009

New benchtop machine offers breakthrough versatility.

Smart, Small, Sleek - A Benchtop Breakthrough

**Filamatic to Demonstrate the New AdaptaFil™ Filling Machine at Interphex
Filamatic Booth #2159**

Baltimore, MD –Filamatic will be demonstrating the AdaptaFil™ Benchtop, a revolutionary semi-automatic filling machine at Interphex March 17th-19th at the Jacob Javits Convention Center in New York. This new Benchtop filler will offer a level of functionality and versatility seen only in fully automatic filling systems and promises to redefine semi-automatic filling operations. The AdaptaFil™ Benchtop's patent pending technology features four easy-changeover metering systems, electronic networking capabilities and sophisticated, user-friendly controls all in a small and sleek stainless steel package about the size of a microwave.

Filamatic will also feature a version of their popular Monobloc fill/finish system. These systems can perform a wide range of tasks such as sorting, feeding, filling, plugging, stoppering, crimping, capping, induction sealing, labeling, and accumulating. The Monobloc is ideal for the Pharma industry because it is engineered to perform multiple tasks within a compact footprint. An optional HEPA filtration system can also be incorporated to provide positive pressure airflow within the machine environment. A variety of liquid metering systems are available to provide a fill volume range of 0.5 ml. to 500 ml. per cycle. In addition, these metering systems are capable of handling a wide range of viscosities. All Filamatic Monoblocs are digitally controlled with menu-driven programs that electronically adjust operating parameters for quick and easy tool-less changeovers. A wide assortment of liquid products can be filled into an array of containers including microtubes, cryovials, glass vials, and plastic containers.

About Filamatic

Filamatic, a division of National Instrument, LLC, manufactures liquid filling equipment and integrated liquid packaging systems. They have been serving the pharmaceutical, biotech, healthcare, cosmetic, chemical specialty and food industries since 1950 and have thousands of liquid filling machines in operation around the world. For more information visit their website at www.filamatic.com or visit the Filamatic booth #2159 at Interphex.

4119-27 Fordleigh Road • Baltimore, Maryland 21215-2292 • Phone (410) 764-0900

<http://www.filamatic.com> • info@filamatic.com • FAX (410) 764-7719